Bruksanvisning Ödemhandskar

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Instruction for use - Edema gloves

Catell recommends that the product is fitted and adapted by a licensed or otherwise qualified health professionals. Users with diabetes, sensitive skin, known contact allergies, impaired sensibility or poor circulation must be extra observant and check for pinching, chafing or edema when using the product. Treatment with compression may be contraindicated in certain conditions e.g. ongoing infection/inflammation, heart problems, tumor disease, kidney and liver disease and peripheral neuropathy. In case of uncertainty consult a physician. Always check material content below. Any serious incidents occurring related to the product should be reported to the manufacturer and to the competent authority of the Member State in which the user is resident. Contact your prescriber or Catell if questions about the product occur.

MODEL	COMPRESSION CLASS	ART.NO	SIZE	SIZE	* DIMENSIONS
Catell Ödemhandske ¾ Finger	1-2 (15-25 mmHg)	5050	XS - L	XXS	13-15 cm
Catell Ödemhandske Fullfinger	1-2 (15-25 mmHg)	6050	XS - L	XS	15-18 cm
EDEMA Light ¾ Finger	1-2 (15-25 mmHg)	7050, 7051	XXS - XL	S	18-20 cm
EDEMA Light Fullfinger	1-2 (15-25 mmHg)	7150, 7151	XXS - XL	M	20-22 cm
EDEMA Light Wide ¾ Finger	1-2 (15-25 mmHg)	7050W	XS - XL	L	22-24 cm
EDEMA Light Long Fullfinger	1-2 (15-25 mmHg)	7160	S - L	XL	24-27 cm
EDEMA Light Thumb	1-2 (15-25 mmHg)	7350	XS - XL		
EDEMA More ¾ Finger	2 (23-32 mmHg)	8050, 8051	XXS - XL		
EDEMA More Fullfinger	2 (23-32 mmHg)	8150, 8151	XXS - XL		

*Measure MCP-joints circumference.
Measurements should be seen as guidelines.
Note - EDEMA is bilateral. Catell Ödemhandske is available in Left/Right.

Properties

- Seams on the outside to minimize pressure on sensitive points.

Use glove

- In purpose to reduce swelling, as recommended by prescriber.
- With the seams outwards.

Contact your prescriber if

- The swelling increases.
- The glove causes numbness.
- The glove causes pain.

Material

Nylon, lycra. Free from latex.

Washing instructions

Follow the washing instructions below. Use a washing pouch. Let dry horizontally.

Single patient use.

Recycling

Sort according to local regulations.

This device conforms to the requirements of the Medical Device Regulation (EU) 2017/745 product class I.

Manufacturer Catell AB

Örkroken 23, SE-138 40 Älta, Sweden | Phone: +46 8 727 95 30 | info@catell.se | www.catell.se 250219:5